



## B and T non-Hodgkin's lymphomas

For B-NHL, Eurofins Biomnis offers five NGS panels:

- the "LLC" NGS panel (see CLL panel data sheet, ref. DS119-INTGB),
- the "LLCTR" NGS panel (see CLL panel data sheet, ref. DS119-INTGB),
- the "LPWAL" NGS panel,
- the "BRAF" NGS panel,
- the "LNHB" NGS panel.

For T-NHL, Eurofins Biomnis offers one panel:

- the "LNHT" NGS panel.

### Non-Hodgkin B lymphoma

- The "LPWAL" NGS panel consists of an analysis of 10 genes: *MYD88*, *CXCR4*, *ARID1A*, *CD79A*, *CD79B*, *NOTCH2*, *TP53*, *BTK*, *PLCG2* and *CARD11*. It must be combined with a bone marrow cytogenetic study.

It is used to diagnose lymphoplasmacytic lymphomas, principally Waldenström's disease. The *MYD88* mutation is reported in more than 90% of cases and can be used to rule out diagnosis of LZM, for example. In more than 40% of cases, a *CXCR4* mutation is observed, generally associated with a *MYD88* mutation. There are also cases with *CXCR4* mutation without *MYD88* mutation. Mutations in the *ARID1A*, *CD79A*, *CD79B* and *NOTCH2* genes have also been reported in Waldenström disease and may contribute to the diagnostic approach.

This panel also provides a prognostic and therapeutic approach: The presence of a *CXCR4* mutation is associated with a poor prognosis and represents a resistance mutation, in the same way as secondary mutations in the *BTK* and *PLCG2* genes.

*TP53* mutation, although rare (2%), is associated with a poor prognosis. The association of a *MYD88* mutation with the absence of a *CXCR4* mutation has a favourable prognosis under BTKi, unlike the association of both *MYD88* and *CXCR4* mutations, which is associated with a late response to BTKi. Absence of *MYD88* mutation in Waldenström disease characterises a disease with a poor prognosis and a poor response to BTKIs or to a combination of bendamustine and rituximab.

The IPSSWM-R prognostic score does not currently include molecular data.

## « LPWAL » NGS Panel - Targeted genes

Gene	Transcript	Exon rank
<i>ARID1A</i>	NM_006015	Full Coding Region
<i>BTK</i>	NM_000061	Full Coding Region
<i>CARD11</i>	NM_032415	Full Coding Region
<i>CD79A</i>	NM_001783	Full Coding Region
<i>CD79B</i>	NM_000626	Full Coding Region
<i>CXCR4</i>	NM_003467	Full Coding Region
<i>MYD88</i>	NM_002468	Full Coding Region
<i>NOTCH2</i>	NM_024408	Full Coding Region
<i>PLCG2</i>	NM_002661	Full Coding Region
<i>TP53</i>	NM_000546	Full Coding Region

**Test code:** LPWAL

- The **BRAF** NGS panel involves mutational analysis of the *BRAF* gene. This mutation is reported in more than 95% of cases of Hairy Cell Leukaemia and represents a very useful diagnostic criterion in association with a targeted immunophenotypic analysis.

## « BRAF » NGS Panel - Targeted genes

Gene	Transcript	Exon rank
<i>BRAF</i>	NM_004333	Full Coding Region

**Test code:** BRAF

- The B-NHL NGS panel consists of an analysis of 45 genes:

*ARID1A/ATM/B2M/BAX/BCL2/BCOR/BIRC3/BRAF/BTK/CARD11/CD79A/CD79B/CREBBP/CXCR4/EGR2/EP300/EZH2/FBXW7/FOXO1A/GNA13/HRAS/KLF2/KRAS/MAP2K1/MCL1/MEF2B/MGA/MYC/MYD88/NFKBIE/NOTCH1/NOTCH2/NRAS/PIM1/PLCG2/POT1/RPS15/SAMHD1/SF3B1/SOCS1/STAT6/TNFAIP3/TNFRSF14/TP53/XPO1*. It must be combined with a blood, bone marrow or lymph node cytogenetic analysis.

Note For CLL, Plasma cell lymphoma and Hairy Cell Leukaemia, see panel sheet DS119-INTGB and corresponding paragraphss.

It is used as a **diagnostic** aid to identify, in addition to histological, cytological and cytogenetic analyses, a type or subtype of B-NHL in the WHO and ICC 2022.

Some examples of **diagnostic** aids:

- For diffuse large cell lymphomas, the GC and ABC subtypes can be characterised by distinct mutational profiles: mutations in *EZH2*, *GNA13*, *MEF2B*, *TNFRSF14*, *B2M* and *CREBBP* genes for the GC type and *MYD88*, *CD79B*, *TNFAIP3*, *CARD11* and *PIM1* genes for the ABC type.
- The CD23+ «follicular lymphoma - non-rearranged BCL2» subtype is generally associated with a *STAT6* or *TNFRSF14* mutation.
- Primary mediastinal large B-cell lymphoma is characterised by mutations in the *STAT6*, *XPO1*, *NFKBIE*, *TNFAIP3*, *GNA13* or *B2M* genes.
- For marginal zone lymphomas, the mutational profile of the NGS panel makes it possible to differentiate extra-nodal forms (*TNFAIP3*, *TNFRSF14* or *TET2* mutations) from nodal or splenic forms (*KLF2* or *NOTCH2* mutations).

This panel also provides **prognostic** support:

- In follicular lymphoma, the m7-FLIPI score is used to assess prognosis using the mutation status of 7 genes (*EZH2*, *ARID1A*, *MEF2B*, *EP300*, *FOXO1*, *CREBBP* and *CARD11*).
- In mantle cell lymphoma, the presence of mutations in *TP53*, *NOTCH1* or *NOTCH2* mutations is associated with a poor prognosis.
- In diffuse large cell lymphoma, new molecular subtypes (MCD, EZB, BN2, ST2, A53 and N1), defined by NGS, can also be used to establish a prognosis.  
This includes the analysis of new genes (*SOCS1*, *NOTCH1* for example). The *TP53* mutation also has a poor prognosis in this entity.

From a **therapeutic** point of view, in addition to research into BTKi resistance mutations (mutations in the *BTK* and *PLCG2* genes) in CLL, MCL, MZL and Hairy Cell Leukaemia, targeted therapies offer new therapeutic choices, for example, the *EZH2* target in follicular lymphoma.

## « B-NHL » NGS Panel - Targeted gene

Gene	Transcript	Exon rank
<i>ARID1A</i>	NM_006015	Full Coding Region
<i>ATM</i>	NM_000051	Full Coding Region
<i>B2M</i>	NM_004048	Full Coding Region
<i>BAX</i>	NM_138761	Full Coding Region
<i>BCL2</i>	NM_000633	Full Coding Region
<i>BCOR</i>	NM_017745	Full Coding Region
<i>BIRC3</i>	NM_001165	Full Coding Region
<i>BRAF</i>	NM_004333	Full Coding Region
<i>BTK</i>	NM_000061	Full Coding Region
<i>CARD11</i>	NM_032415	Full Coding Region
<i>CD79A</i>	NM_001783	Full Coding Region
<i>CD79B</i>	NM_000626	Full Coding Region
<i>CREBBP</i>	NM_004380	Full Coding Region

Gene	Transcript	Exon rank
<i>CXCR4</i>	NM_003467	Full Coding Region
<i>EGR2</i>	NM_000399	Full Coding Region
<i>EP300</i>	NM_001429	Full Coding Region
<i>EZH2</i>	NM_004456	Full Coding Region
<i>FBXW7</i>	NM_033632	Full Coding Region
<i>FOXO1A</i>	NM_002015	Full Coding Region
<i>GNA13</i>	NM_006572	Full Coding Region
<i>HRAS</i>	NM_176795	Full Coding Region
<i>KLF2</i>	NM_016270	Full Coding Region
<i>KRAS</i>	NM_033360	Full Coding Region
<i>MAP2K1</i>	NM_002755	Full Coding Region
<i>MCL1</i>	NM_021960	Full Coding Region
<i>MEF2B</i>	NM_001145785	Full Coding Region

Gene	Transcript	Exon rank
<i>MGA</i>	NM_001164273	Full Coding Region
<i>MYC</i>	NM_002467	Full Coding Region
<i>MYD88</i>	NM_002468	Full Coding Region
<i>NFKBIE</i>	NM_004556	Full Coding Region
<i>NOTCH1</i>	NM_017617	Full Coding Region
<i>NOTCH2</i>	NM_024408	Full Coding Region
<i>NRAS</i>	NM_002524	Full Coding Region
<i>PIM1</i>	NM_002648	Full Coding Region
<i>PLCG2</i>	NM_002661	Full Coding Region
<i>POT1</i>	NM_015450	Full Coding Region

Gene	Transcript	Exon rank
<i>RPS15</i>	NM_001018	Full Coding Region
<i>SAMHD1</i>	NM_015474	Full Coding Region
<i>SF3B1</i>	NM_012433	Full Coding Region
<i>SOCS1</i>	NM_003745	Full Coding Region
<i>STAT6</i>	NM_003153	Full Coding Region
<i>TNFAIP3</i>	NM_006290	Full Coding Region
<i>TNFRSF14</i>	NM_003820	Full Coding Region
<i>TP53</i>	NM_000546	Full Coding Region
<i>XPO1</i>	NM_003400	Full Coding Region

**Test code:** LNHB

## Non-Hodgkin T lymphoma

- The "LNHT" NGS panel consists of an analysis of 19 genes: *ARID1A*, *ATM*, *BCOR*, *CARD11*, *CD28*, *DNMT3A*, *EP300*, *FBXW7*, *IDH2*, *JAK2*, *JAK3*, *MGA*, *NOTCH1*, *PLCG1*, *RHOA*, *STAT3*, *STAT5B*, *TET2* and *TP53*. It must be combined with a search for T clonality and a cytogenetic study.

It is primarily intended as a **diagnostic** aid for peripheral T lymphomas. Recurrent cytogenetic abnormalities are also reported for diagnostic purposes and must be investigated in association with these abnormalities. For example, the combination of *IDH2*, *RHOA* and *TET2* mutations confirms the diagnosis of angioimmunoblastic TFH lymphoma (the *IDH2* mutation being exclusive to this subtype of TFH lymphoma). The presence of a *STAT5B* and *JAK3* mutation helps in the differential diagnosis between a monomorphic epitheliotropic intestinal T-cell lymphoma (MEITL) (*STAT5B* mutation reported in more than 60% of cases) and an enteropathy-associated T-cell lymphoma (EATL), which tends to be associated with a *STAT3* mutation (20% of cases). Conversely, ALK+ anaplastic large cell lymphoma is associated with virtually no mutations in the classic T-cell lymphoma mutational spectrum.

In terms of **theranostics**, potential targeted therapies are currently being studied, targeting the *IDH2* or *RHOA* genes, for example.

## Gene mutation frequencies and mutations of high diagnostic interest in T-cell lymphomas

	T-PLL	T-LGLL	ATLL	SS	EATL	MEITL	HSTL	ALCL	TFHL-AI	ENTKL
ARID1A				X						X
ATM	X			X						
BCOR										X
CARD11			X	X					X	
CD28			X	X					X	
DNMT3A				X					X	
EP300				X						X
FBXW7			X							
IDH2									X	
JAK2										X
JAK3	X					X		X		X
MGA										X
NOTCH1			X							
PLCG1			X	X					X	
RHOA									X	
STAT3		X	X	X	X		X	X(ALK-)		X
STAT5B	X	X		X	X	X	X			X
TET2	X				X				X	X
TP53	X		X	X	X	X		X		X

### Frequency of mutation:

X 5% < < 20%

X 20% < < 50%

X 50% < < 100%

■ Mutation of high diagnostic interest

### Abbreviations:

T-PLL: T-prolymphocytic leukaemia

T-LGLL: T-large granular lymphocytic leukaemia

ATLL: Adult T-cell leukaemia/lymphoma

SS: Sezary syndrome

EATL: Enteropathy-associated T-cell lymphoma

MEITL: Monomorphic epitheliotropic intestinal T-cell lymphoma

HSTL: Hepatosplenic T-cell lymphoma

ALCL: Anaplastic large cell lymphoma

TFHL-AI: T-follicular helper (TFH) cell lymphoma, angioimmunoblastic-type

ENTKL: Extranodal NK/T-cell lymphoma

## « LNHT » NGS Panel - Targeted gene

Gene	Transcript	Exon rank
<i>ARID1A</i>	NM_006015	Full Coding Region
<i>ATM</i>	NM_000051	Full Coding Region
<i>BCOR</i>	NM_017745	Full Coding Region
<i>CARD11</i>	NM_032415	Full Coding Region
<i>CD28</i>	NM_006139	Full Coding Region
<i>DNMT3A</i>	NM_022552	Full Coding Region
<i>EP300</i>	NM_001429	Full Coding Region
<i>FBXW7</i>	NM_033632	Full Coding Region
<i>IDH2</i>	NM_002168	Full Coding Region
<i>JAK2</i>	NM_004972	Full Coding Region
<i>JAK3</i>	NM_000215	Full Coding Region
<i>MGA</i>	NM_001164273	Full Coding Region
<i>NOTCH1</i>	NM_017617	Full Coding Region
<i>PLCG1</i>	NM_002660	Full Coding Region
<i>RHOA</i>	NM_001664	Full Coding Region
<i>STAT3</i>	NM_139276	Full Coding Region
<i>STAT5B</i>	NM_012448	Full Coding Region
<i>TET2</i>	NM_001127208	Full Coding Region
<i>TP53</i>	NM_000546	Full Coding Region

**Test code:** LNHT

**Pre-analytical requirements:** EDTA whole blood or bone marrow

**Turnaround time :** 13 days (Results may require an extended turnaround time of one week, depending on the confirmation tests required by Sanger sequencing).

## Contact

Eurofins Biomnis Ireland  
sales@ctie.eurofinseu.com  
Phone: 1800 303 349

## References

The 5th edition of the World Health Organization Classification of Haematolymphoid Tumours: Lymphoid Neoplasms. Alaggio R et al. Leukemia. 2022 Jul;36(7):1720–1748. PMID: 35732829

The International Consensus Classification of Mature Lymphoid Neoplasms: a report from the Clinical Advisory Committee. Campo E et al. Blood. 2022 Sep 15;140(11):1229–1253. PMID: 35653592

Diagnostic and prognostic molecular pathology of lymphoid malignancies. Fend F et al. Virchows Arch 2023 Sep 25. PMID: 37747559

Genomic profiling for clinical decision making in lymphoid neoplasms. De Leval L. et al. Blood. 2022 Nov 24;140(21):2193–2227. PMID: 36001803

Precision diagnostics in lymphomas – Recent developments and future directions. Mansouri L et al. Semin Cancer Biol. 2022 Sep;84:170–183. PMID: 34699973

Genetic Landscape of Peripheral T-Cell Lymphoma. Hathuc V, Kreisel F. Life (Basel). 2022 Mar 11;12(3):410. PMID: 35330161

Molecular profiling in diffuse large B-cell lymphoma: why so many types of subtypes. Morin R et al. Br J Haematol. 2022 Feb;196(4):814–829. PMID: 34467527

Exploiting gene mutations and biomarkers to guide treatment recommendations in mantle cell lymphoma. Goy A. Expert Rev Hematol. 2021 Oct;14(10):927–943. PMID: 34253131

Genomic Landscape of Waldenström Macroglobulinemia and Its Impact on Treatment Strategies. Treon S. et al. J Clin Oncol. 2020 Apr 10;38(11):1198–1208. PMID: 32083995

A revised international prognostic score system for Waldenström's macroglobulinemia. Kastiris et al. Leukemia 2019 Nov;33(11):2654–2661. PMID: 31118465

The Need for a Consensus Next-generation Sequencing Panel for Mature Lymphoid Malignancies. Sujobert P et al. Hemasphere. 2018 Dec 27;3(1):e169. PMID: 31723808

Integration of gene mutations in risk prognostication for patients receiving first-line immunochemotherapy for follicular lymphoma: a retrospective analysis of a prospective clinical trial and validation in a population-based registry. Pastore A et al. Lancet Oncol. 2015 Sep;16(9):1111–1122. PMID: 26256760