

## **TUH - Sample Packing Guide**

Eurofins Biomnis Ireland (Clinical Chem Samples)

## **Urgent Samples**

Red Bag

Potassium		Phosphate		
Potassium needs to get to us within 2 hrs of venepuncture (unless spun)		Phosphate needs to get to us within 6 hrs of venepuncture (unless spun)		
Please arrange phlebotomy times in line with your Logistics collection times.				
•	For Urgent Requests, please provide clinicial reasons along with Out of Hours contact number and name of the Requesting Clinician. In the absence of these details, the test request will be processed as routine.			
CONTACT POINT	email: clientservices@ctie.eurofinseu.com for an issues regarding results for the Eurofins Biomnis tests.			



Samples for Eurofins Biomnis Ireland are to be packed separately in the Blue/Red bags addressed to EBI.



## **Tallaght University Hospital (TUH) Samples**

Samples for TUH are denoted in the "Laboratory" section of the form, depending on their department.				
Laboratory: TUHBio		Laboratory: TUHHaem		
Laboratory: TUHImm		Laboratory: TUHMicro		
•	These are to be packed separately to Eurofins Biomnis Ireland samples, as per usual TUH procedures.			
CONTACT POINT	For detailed information on TUH requirements please refer to <a href="https://tuh.ie/Healthcare-Professionals/Contacts-Information-/Laboratory-Medicine.html">https://tuh.ie/Healthcare-Professionals/Contacts-Information-/Laboratory-Medicine.html</a> Email: <a href="mailto:gplabqueries@tuh.ie">gplabqueries@tuh.ie</a> for an issues regarding results for the TUH tests.			

MK05 Issue No. 1.01 Active Date: 31/10/25